

MARPE

Healthy Life Leaders


GTS (Gait Training System)

- Smooth weight assisting rail system: Walking rail, gait equipment, suit, and others
- Enables independent and stable gait training of users
- Enables efficient exercise for the activities of daily living (ADL) of senior citizens, rehabilitation patients or person with disabilities in hospitals, rehabilitation center, and other facilities
- A group therapeutic exercise that allows multiple users to use it at the same time
- Enables step-by-step circuit rehabilitation training indoors for adaptation to walking outdoor
- Medical device manufacturing report: Manual orthopedic exercise device (MGTS-00101)

M-Sling

It is a movable sling that can change direction and move 360 degrees. In addition, various fixed sling exercise programs are possible by fixing them to the rail. So both walking and fixed sling exercises can be performed, which is useful for economic benefits and space utilization.


357mm(Width) x 65mm(Depth) x 120mm(Height) / 5kg



Single-Sling

It is a product that can be moved along the rail and can be fixed to the rail for limited fixed sling exercise.

130mm(Width) x 65mm(Depth) x 165mm(Height) / 2.1kg



Manual


Auto

TTS

"Turn Terminal" is used to change the direction of the sling when walking as well as parking, and electric models can rotate not only in T-shape but also in a + - shape.

993mm(Width) x 696mm(Depth) x 110mm(Height) / 21kg

GTS-M Sling



GTS-Single Sling



TTS (Turn Terminal System)



GTS-Ceiling type


“GTS-Ceiling type” is a ceiling-mounted product that enables customized rehabilitation training and gait training by providing various shapes and sizes at the request of clients.



GTS-Stand type

“GTS-Stand type” is an assembly type product that can be installed without any on-site construction, which enables various rehabilitation training and gait training using parallel bars.

4000mm(Width) x 2500mm(Depth) x 2400mm(Height) / 300kg



GTS-Basic



GTS-Advanced



GTS-Circuit



GTS-Stand

