Gait Training System-Rail

GTS-Ceiling Type



"GTS-Ceiling Type" is a ceiling-mounted product that enables customized rehabilitation training and gait training by providing various shapes and sizes at the request of clients.

GTS-Stand Type



"GTS-Stand Type" is an assembly-type product that can be installed without any on-site construction, which enables various rehabilitation training and gait training using parallel bars.


Features of Rail

- The rail design prevents the Sling from slipping unnecessarily, enabling smooth gait rehabilitation.
- It can be customized to the size and shape desired by the consumer.
- No problem with long-term use because it is a safe and robust installation.

Specifications

Customized in various shapes.

Details



Customizing

We offer custom rails in the desired size and shape.



Rail bending

The only company in Korea that does rail bending directly, so the quality is excellent.



CE Certification

Certified with European CE certification for consumer safety and health.



Parallel bar

The stand-type product includes a Parallel Bar for professional gait rehabilitation.









Gait Training System-Turn Terminal

- "Turn Terminal System (TTS)" is used to change the direction of the sling when walking as well as it is used for parking the unused sling device
- In the case of auto model, T-shaped rotation as well as +-shaped rotation is possible

TTS-Manual



TTS-Auto



Features of Turn Terminal

- It is used for direction change and sling storage when configuring the GTS rail shape.
- Can choose between manual and auto type, and both models provide smooth direction change.

Specifications

- Width 1000mm
- Height 110mm
- Depth 700mm
- Weight 25kg

Details



Manual handle

Rotate with one handle for convenient operation



Motorized device

Can change the direction comfortably with the remote control, reducing user fatigue.


CE Certification

Certified with European CE certification for consumer safety and health.



Customizing

For efficient use, we provide suitable TTS for your facility.









Gait Training System-M Sling

- M Sling is a movable sling that can change direction and rotate 360 degrees
- Various sling exercise programs are possible by fixing to the rail



Features of M Sling

- As a representative sling product of GTS, it is the most functional product.
- Using the principle of the sling, not only walking training but also upper and lower body strength training, relaxation training, stretching, and balance training are possible.

Specifications

- Width 357mm
- Height 120mm
- Depth 65mm
- Weight 5kg

Details



Roller

The specially designed roller device enables smooth movement and free change of direction.



Fixing device

It is used to securely fixed the M Sling in the desired location.

CE Certification

Certified with European CE certification for consumer safety and health.



M Suit

It is easy to use, and It is a suit made of only high-quality material for the comfort of the user.









Gait Training System-Single Sling

- Single Sling is a movable sling that allows free rail movement
- A balance training program is possible by fixing it to the rail



Features of Single Sling

- Single sling is specialized for gait training.
- The size and weight are reduced to enable a smoother gait during gait training.

Specifications

- Width 130mm
- Height 165mm
- Depth 65mm
- Weight 2.1kg

Details



Roller

The specially designed roller device enables smooth movement and free change of direction.



Fixing device

It is used to securely fixed the Single Sling in the desired location.

CE Certification

Certified with European CE certification for consumer safety and health.



Accessories

We provide the customized holders and accessories for various rehabilitation training.







